

Histological Slide Preparation

By:

Dr. Najlaa Al-Mansour

preparation of tissue for histological slide:

- ❖ Fixing tissue.
- ❖ Embedding tissue.
- ❖ Sectioning the specimen.
- ❖ Basic Staining technique.

1-Fixing tissue:

- ❖ The specimens are delivered to the Histology Laboratory immediately after the autopsy.
- ❖ The tissue is labeled with an identifying Medical Examiner's case number and is placed in plastic capsules or cassettes.
- ❖ The fixative used to preserve tissue and maintenance life like structure and prevent decomposition .



2- Embedding

- ❖ **Preparation:** Cut (1Cm.) of specimen and put in cassette, then the cassette is placed in a 10% **formalin** fixative solution,(24-48 hours).
- ❖ **Tissue processing:** dehydration, clearing, infiltration and blocking.
 - dehydrate the specimen with ascending series of alcohol(70%-90%-100%).
 - clearing with organic solvent e.g. (xylene)which removes alcohol then allowed infiltration by paraffin wax.
 - Length of processing depends on size and density of specimen (hours to overnight).
- ❖ **Blocking:** remove the cassette and put it on warm plate, then the specimen removed with warm forceps and placed in the mould and filled with wax.

3- Sectioning the specimen:



- ❖ Using rotary microtome to cut the specimen into (3-10micron).
- ❖ Float ribbon (tissue)on warm water to flatten then lifts the paraffin ribbon on the microscopic glass slide .
- ❖ Allow to dry at 37C.to ensure section to adhere to slide.



4-Staining:


- ❖ Rehydration: descending series of alcohol (100%-90%-70%).
- ❖ Rinse with tap water.
- ❖ Staining with basic stain e.g. hematoxylin and eosin, the hematoxylin staining the nuclei blue (10 min.), then the slide wash with tap water to neutralized the ph.
- ❖ Using eosin to stain the cytoplasm, collagen and muscle fiber.
- ❖ Dehydrate the slide with alcohol and use the xylene (15min.) to removed the traces of alcohol.
- ❖ Mounting the slide by using (DPX) then put the cover slip.

Basic Stain



Special Staining

- ❖ Special staining can be used for specific purposes.
- ❖ Special stains are used to determine bacteria, virus, fungi, intercellular, and intracellular structures.
- ❖ E.g. Masson Trichrom stain and Silver stain.



Signing Out

- The last step includes entering the case into the computer database and printing the labels.

